

DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES

Enhancing intracellular delivery of a **broad range** of proteins

USER PROTOCOL – #DIV031F1

ABOUT THE NANOPARTICLES	1
OVERVIEW	1
COMPONENTS	1
STORAGE	1
EQUIPMENT AND MATERIALS REQUIRED BUT NOT SUPPLIED	2
CONSIDERATIONS BEFORE STARTING	2
DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES	3
PROTEIN MODIFICATION STEP	3
PROTEIN PURIFICATION STEP	4
DIVERSA FORMULATION STEP	6
DIVERSA FLUOGREEN PROTEIN FORMULATION STEP	7
EXAMPLE OF β -GAL ASSOCIATION PROTOCOL.....	9
EXAMPLE OF UPTAKE ASSAY PROTOCOL	10
OPTIMIZATION GUIDELINES.....	11
RECOMMENDATIONS OF USE AND TECHNICAL NOTES	12
FREQUENTLY ASKED QUESTIONS	14
ONLINE RESOURCES	14

ABOUT THE NANOPARTICLES

OVERVIEW

DIVERSA is a biocompatible and biodegradable cell-friendly technology for enhancing intracellular delivery of proteins, paving the way towards clinical translation.

DIVERSA PROTEIN DELIVERY NANOPARTICLES, based on lipids carrying a reactive group, is suitable for click chemistry. A modification in your protein/s of interest is needed for the association to the **DIVERSA NANOPARTICLES** without comprising its structure and activity. Azido (N₃- group)-modified proteins will covalently react with the group exposed on the surface of the **DIVERSA NANOPARTICLES** by following easy and mild conditions.

DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES uses strongly labelled fluorescent nanoemulsions that are easily internalized by live cells and can be visualized by a wide variety of fluorescent platforms (flow cytometry, microplate assays, fluorescence, and confocal microscopy) in less than two hours at Ex/Em = 495/503 nm.

COMPONENTS

- 1x **DIV-LINKER**.
- 1x **DIV031F1** vial for reconstitution.
- 1x **DIVTECH** vial for preparation of **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES**.
- 1x Ultracentrifugation filter 0.5 mL- 10 kDa.
- 2x Colector microtubes.
- 2x Tips for 1 mL micropipette.

STORAGE

Before formulating, store the vials at -20 °C. Once formulated, storage is recommended at 4 °C.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EQUIPMENT AND MATERIALS REQUIRED BUT NOT SUPPLIED

- 1 mL micropipette.
- Sterile standard microtubes.
- MilliQ water or any other recommended buffer.
- Ethanol (EtOH) 96%.
- PBS 2 mM at pH 7.4.
- DMSO.
- Protein/s of interest.

CONSIDERATIONS BEFORE STARTING

- The proteins need to be functionalized with the **DIV-LINKER** for an efficient association to **DIVERSA FLUOGREEN NANOPARTICLES** (follow the protocol below).
- The following protocol is optimized for the preparation of 1.5 mL of **DIVERSA FLUOGRENN PROTEIN DELIVERY NANOPARTICLES** (starting from one **DIV031F1** vial for reconstitution).
- **DIVERSA** cannot guarantee the optimal characteristics of the final formulation if modifications in the protocol are introduced.
- It is recommended to use **DIVERSA FLUOGREEN NANOPARTICLES** within 60 days.
- **DIVERSA FLUOGREEN NANOPARTICLES** is stable with supplemented cell culture media at 37 °C: DMEM, RPMI.
- Do NOT use any buffer solution containing Triton-X, SDS or Tween-20 for the preparation or manipulation of **DIVERSA/DIVERSA FLUOGRENN PROTEIN DELIVERY NANOPARTICLES**.
- Do NOT freeze **DIVERSA/DIVERSA FLUOGRENN PROTEIN DELIVERY NANOPARTICLES**.
- Do NOT heat up **DIVERSA/DIVERSA FLUOGRENN PROTEIN DELIVERY NANOPARTICLES**.
- Do NOT centrifuge or vortex **DIVERSA/DIVERSA FLUOGRENN PROTEIN DELIVERY NANOPARTICLES**.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES PROTOCOL

PROTEIN MODIFICATION STEP:

Note: pH 7.2 - 8 is recommended for protein modification.

1. Add 40 μ L of DMSO to the **DIV-LINKER** microtube. Vortex gently and then, spin-down.
2. Select the volume of **DIV-LINKER** microtube based on the MW and the amount of your protein as we recommend in [Table 1](#) (Recommendations of Use and Technical Notes). Then, dilute the **DIV-LINKER** solution up to 250 μ L in PBS 2 mM at pH 7.4 or the buffer suggested in [Table 2](#) (Recommendations of Use and Technical Notes).

Note: Avoid buffers containing b-mercaptoethanol, ammonium salts and primary amines.

3. Add the protein of interest to the **DIV-LINKER** microtube and adjust the final volume to 350 μ L with PBS 2 mM at pH 7.4 or the buffer suggested in [Table 2](#) (Recommendations of Use and Technical Notes). Vortex gently and then, spin-down.
4. Incubate the **PROTEIN-LINKER** for 6 h at room temperature (RT) and then, incubate overnight at 4 $^{\circ}$ C.

1. Add 40 μ L of DMSO to the **DIV-LINKER** microtube
2. Select the volume of **DIV-LINKER** vial and dilute it up to 250 μ L
3. Add the protein to the **DIV-LINKER** vial. Adjust the volume to 350 μ L
4. Incubate at RT for 6 h, and at 4 $^{\circ}$ C overnight




Figure 1. **PROTEIN MODIFICATION STEP** protocol.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

PROTEIN PURIFICATION STEP:

5. Equilibrate the Amicon Ultracentrifugation filter 0.5 mL- 10 kDa with PBS 2 mM at pH 7.4 or the buffer used in the previous protocol. Centrifuge at 14,000 RCF for 5 min at 4 °C. Discarding the flow-through and the filter is ready to use.
6. After incubation time, transfer the **PROTEIN-LINKER** to the Amicon Ultracentrifugation filter 0.5 mL- 10 kDa, and top up the volume to 500 µL with PBS 2 mM at pH 7.4 or the buffer suggested in [Table 2](#) (Recommendations of Use and Technical Notes).
7. Centrifuge at 14,000 RCF for 5 min at 4 °C.

Note: The volume should drop to ~ 300 µL.

8. Discard the flow-through. Wash the **PROTEIN-LINKER** four times with PBS 2 mM at pH 7.4 or the buffer suggested in [Table 2](#) (Recommendations of Use and Technical Notes) as performed in the previous step. Centrifuge at 14,000 RCF for 5 min at 4 °C.

Note: the centrifugation time of the last wash depends on the desired volume you want to concentrate your protein.

Note: we recommend concentrate the **PROTEIN-LINKER** at 2.5 mg/mL.

9. Invert the filter containing the **PROTEIN-LINKER** into a clean collector microtube, and centrifuge at 1,000 RCF for 2 min at 4 °C. A 100% of **N3-MODIFIED PROTEIN** recovery is guaranteed by this method.

The **N3-MODIFIED PROTEIN** is now ready-to-use. Alternatively, keep it at -20 °C and use it in the following 3 months.


Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

6. Transfer the **PROTEIN-LINKER** to the pre-equilibrated filter (provided).
7. Centrifuge at 14,000 RCF for 5 min at 4 °C
8. Discard the flow-through and wash four times as performed previously.



9. Invert the filter containing the **PROTEIN-LINKER** into a clean collector microtube. Centrifuge at 1,000 RCF for 2 min at 4 °C.




Figure 2. PROTEIN PURIFICATION STEP protocol.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

DIVERSA FORMULATION STEP:

10. Reconstitute the **DIV031F1** vial with 100 μL of EtOH. Pipette up and down gently for mixing the lipids trying to recover all of them from the wall vial and keep the suspension in the vial.
11. Add 900 μL of ultrapure water into the **DIVTECH** vial.
12. Transfer the whole volume from **DIV031F1** vial to the **DIVTECH** vial, using a micropipette and the 1 mL micropipette tip provided.

IMPORTANT: Before adding the lipids from **DIV031F1** vial to **DIVTECH** vial, set the micropipette at the maximum volume to have dead air volume in the tip for mixing in a faster and vigorous way. Then, place the 1 mL micropipette tip into the buffer solution of **DIVTECH** vial, and pipette up and down for 30 seconds, avoiding any spillage.

The **DIVERSA FLUOGREEN NANOPARTICLES** is now ready for the association of the protein. Alternatively, keep it at 4 °C and use it in the following 60 days.




Figure 3. DIVERSA FLUOGREEN FORMULATION STEP protocol.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

DIVERSA FLUOGREEN PROTEIN ASSOCIATION STEP:

13. See the following [Table](#) to prepare **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** depending on the required quantity of your protein to obtain a maximum yield of association of 6.5 μM :

DIVERSA FLUOGREEN NANOPARTICLES	N3-PROTEIN solution	μg of $\beta\text{-GAL}$
1 mL	500 μL	1755 μg
500 μL	250 μL	877.5 μg
200 μL	100 μL	351 μg
100 μL	50 μL	175.5 μg
50 μL	25 μL	87.75 μg

Note: Follow this equation to calculate the required quantity of your protein with the recommended concentration at 6.5 μM .

$$\text{Protein } (\mu\text{g}) = 6.5 \mu\text{M} \times \mathbf{V} (\mathbf{N3} - \mathbf{PROTEIN} \text{ solution}) \times \text{Mw} (\text{kDa}) \times 10^{-3}$$

Note: For example, to associate 50 μL of N3-modified $\beta\text{-Gal}$ (540 kDa) with 100 μL of **DIVERSA FLUOGREEN NANOPARTICLES** to obtain a maximum yield of protein associated formulation (6.5 μM), you should add 175.5 μg of N3-modified $\beta\text{-Gal}$.

$$\text{N3} - \beta - \text{Gal} = 6.5 \mu\text{M} \times 50 \mu\text{L} \times 540 \text{ kDa} \times 10^{-3} = 175.5 \mu\text{g}$$

14. Add the **N3-MODIFIED PROTEIN** into **DIVERSA FLUOGREEN NANOPARTICLES**. Pipette up and down gently for 30 seconds, avoiding any spillage.
15. Incubate the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** for 4 h at RT and then, incubate overnight at 4 $^{\circ}\text{C}$.

The **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** is now ready-to-use. Alternatively, keep it at 4 $^{\circ}\text{C}$ and use it in the following 30 days.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.




Figure 1. **DIVERSA FLUOGREEN PROTEIN FORMULATION STEP** protocol.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EXAMPLE OF β -GAL ASSOCIATION PROTOCOL

1. Add 40 μ L of DMSO to the **DIV-LINKER** microtube. Vortex gently and then, spin-down.
2. Following the [Table 1](#), to associate 175 μ g of β -Gal (540 kDa), take 5 μ L of the **DIV-LINKER** microtube and dilute up to 250 μ L in dH₂O.
3. Add 35 μ L of β -Gal solution (stock concentration at 5 mg/mL) to the **DIV-LINKER** microtube and adjust the final volume to 350 μ L with dH₂O. Vortex gently and then, spin-down.
4. Incubate the **PROTEIN-LINKER** for 6 h at RT and then, overnight at 4 °C.
5. After incubation time, transfer the **PROTEIN-LINKER** to the pre-equilibrated filter and add dH₂O to up 500 μ L. Centrifuge at 14,000 RCF for 5 min at 4 °C.
6. Discard the flow-through and wash the **PROTEIN-LINKER** four times with dH₂O centrifuging at 14,000 RCF for 5 min at 4 °C.
7. Invert the filter into a clean collector microtube and centrifuge at 1,000 RCF for 2 min at 4 °C to collect the 100% of the β -Gal-N3.
8. Reconstitute the **DIV031F1** vial with 100 μ L of EtOH. Pipette up and down gently for mixing the lipids.
9. Add 900 μ L of ultrapure water into the **DIVTECH** vial.
10. Add the content of **DIV031F1** to the **DIVTECH** vial.
11. Transfer the whole volume from **DIV031F1** vial to the **DIVTECH** vial, using a micropipette and the 1 mL micropipette tip provided. The **DIVERSA FLUOGREEN NANOPARTICLES** is now ready for the association of the protein.
12. To associate 175 μ g of β -Gal-N3 (540 kDa) and following the previous [Table](#): add 50 μ L of the β -Gal-N3 into 100 μ L of the **DIVERSA FLUOGREEN NANOPARTICLES**. Pipette up and down gently for 30 seconds, avoiding any spillage.
13. Incubate the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** for 4 h at RT followed by overnight at 4 °C.

The **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** is ready-to-use or alternatively, keep it at 4 °C and use it in the following 30 days.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EXAMPLE OF UPTAKE ASSAY PROTOCOL

1. Seed the recommended number of the cells in a 6-well plate the day before of the uptake experiment for FACS analysis.

Note: the cell density should be determined for each cell type and well plate.

2. Prepare the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** following the provided protocol.

3. Add 960 μ L of fresh cell cultured medium supplemented with 10% (v/v) of FBS and, if necessary, 1 % (v/v) of antibiotics.

4. Add the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** to a final volume of 1 mL.

Note: The final concentration of the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** is calculated considering the final volume of 1 mL.

5. Incubate the cells for 2-4 hours at 37 °C.

Note: Depending on the type of readout assay performed, incubation times may influence delivery efficiency.

7. After incubation time, remove the cell culture medium with the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** and, carefully, wash the cells twice with PBS 1X buffer and remove it.

The **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** is efficiently internalized inside the cells.

Note: we recommend wash the cells with PBS 1X buffer containing calcium and magnesium ions to avoid maximum detachment of living cells.

8. Adjust the concentration for the analysis using the flow cytometer and the following parameters: wavelength of excitation at 595 nm, wavelength of emission at 503 nm.

Note: we recommend analysis on the same day of the experiment. However, for extended storage (> 16 h), we recommend resuspend the cells in 4% (v/v) of paraformaldehyde to prevent cell deterioration.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

OPTIMIZATION GUIDELINES

It is highly recommended to optimize your conditions to get the best **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** performance. Optimize one parameter at a time.

The following parameters can be optimized:

- **Amount of DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES:** Start fixing the concentration and amount of your protein to be delivered, and then you may vary the quantity of the **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** used.
- **Amount of protein to be delivered:** you may need to vary the amount of your protein to be delivered. Depending on the sensitivity of your assay, a greater amount of protein and **DIVERSA FLUOGREEN PROTEIN DELIVERY NANOPARTICLES** may be required.
- **Concentration of the protein solution:** We recommend protein concentration at 6.5 μM . At lower concentrations, we recommend concentrating your protein using ultracentrifugation filters 0.5 mL- 10 kDa. Use the corresponding buffer in [Table 2](#) (Recommendations of Use and Technical Notes).
If your protein is in powder, we recommend dissolve it at final concentration of 6.5 μM , using the corresponding buffer in [Table 2](#) (Recommendations of Use and Technical Notes).
- **Buffer used in PROTEIN PURIFICATION STEP:** you can try different buffers as suggested in [Table 2](#) (Recommendations of Use and Technical Notes). For some biomolecules, the buffer used may be critical, for example PBS buffer works well with β -Galactosidase, but not the Tris buffer.
- **Cell type and density:** you may need to optimize cell numbers. Delivery efficacy may be sensitive to the confluency of the cells in culture.
- **Incubation times for *in vitro* assays:** you may vary incubation times, depending on the type of functional assay performed, shorter or longer incubation time may influence delivery efficiency.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

RECOMMENDATIONS OF USE AND TECHNICAL NOTES

Table 1. Suggested volume of the DIV-LINKER solution based on the required amount and M_w of the protein of interest.

M_w PROTEIN	PROTEIN	DIV-LINKER	BUFFER
> 250 kDa	> 50 μ g	5 μ L	Up to 350 μ L
	20-50 μ g	2 μ L	
	1-20 μ g	1 μ L	
15-250 kDa	> 50 μ g	40 μ L	Up to 350 μ L
	20-50 μ g	15 μ L	
	1-20 μ g	5 μ L	
< 15 kDa	> 50 μ g	40 μ L	Up to 350 μ L
	20-50 μ g	20 μ L	
	1-20 μ g	10 μ L	

Table 2. Suggested buffer for the PROTEIN MODIFICATION STEP and PROTEIN PURIFICATION STEP for *in vivo* experiments.

BUFFER	CONCENTRATION
Physiological saline	0.9% NaCl
Glucose solution for infusion	5-10% w/v
Sodium DL-Lactate solution	60% w/w

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

Table 3. Suggested buffer solutions for buffer exchange after **DIVERSA FORMULATION**.

BUFFER SOLUTION	CONCENTRATION
Ultrapure water	N/A
PBS	2-50 mM
NaCl	150 mM
HEPES	10-25 mM
MES	10-50 mM
Tris-HCl	0.1-1M

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

FREQUENTLY ASKED QUESTIONS

QUESTION	ANSWER
What is the role of the modification of the protein?	The modification of the proteins facilitates their association to DIVERSA FLUOGREEN NANOPARTICLES by following easy and mild procedures that do not compromise its structure and activity. DIVERSA does not include cationic components.
Can the protein modification alter its biological function?	Usually, the proposed modifications do not alter the native protein biological function. Please contact DIVERSA for advice depending on your specific molecule. DIVERSA can provide alternative procedures.
Are there additional modifications that can be suitable to my protein?	Please contact DIVERSA for advice depending on your specific molecule.
How can I measure the size of the final formulation?	Diameter size can be measured by Dynamic Light Scattering (DLS) analysis adding to the cuvette 20 µL of DIVERSA FLUOGREEN PROTEIN NANOPARTICLES with 180 µL of milliQ water.
What if I need to work with higher protein concentrations than the ones provided in Table 2?	You can concentrate the formulation (see next question), or alternatively, contact DIVERSA for advice depending on your specific protein.
How do I concentrate the formulation?	If necessary, the final volume of DIVERSA/DIVERSA FLUOGREEN PROTEIN NANOPARTICLES can be concentrated by using a SpeedVac or Rotavap in mild conditions (avoid surpassing 35 °C or drying out the samples). Samples can be concentrated up to 4-fold its original volume (i.e., to a final volume 250 µL).

ONLINE RESOURCES

Visit our website www.diversatechnologies.com for further information.

Shipping temperature may differ from storage temperature. This does not alter the performance of the product.

DIVERSA TECHNOLOGIES S.L. | Edificio Emprendia, Campus Sur, 15782, Santiago de Compostela, Spain.

Technical support: email: info@diversatechnologies.com | www.diversatechnologies.com

FOR RESEARCH USE ONLY. Not for use in humans. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.